


BenchTop Plate Pipettor

Independent Test Results

Well Volume Goal: 10 μ l

	1	2	3	4	5	6	7	8	9	10	11	12
A	10.06	10.08	10.03	10.06	10.01	9.99	10.03	9.98	10.01	9.87	9.93	9.83
B	10.04	10.09	10.01	10.20	10.09	10.14	9.99	10.03	10.01	9.92	9.86	9.79
C	9.94	10.08	10.01	9.96	9.94	9.91	9.95	9.91	10.01	9.97	9.95	9.90
D	9.99	9.97	9.97	9.97	10.01	9.91	9.94	9.84	10.01	9.96	9.95	9.96
E	10.02	10.08	9.99	9.98	9.91	9.91	9.99	10.09	10.01	10.03	10.02	9.89
F	10.19	10.06	10.14	9.95	9.90	9.97	9.93	10.03	10.01	10.09	10.12	9.97
G	10.01	10.15	10.07	9.99	9.89	9.97	10.02	9.88	10.01	9.98	10.08	9.95
H	10.17	10.18	10.10	10.03	10.07	10.10	10.17	10.10	10.01	10.02	9.99	10.06



Overall Results: Volume 10 μ l

Avg. Volume	10.01 μ l
Standard Deviation	0.09 μ l
Accuracy	0.06% Pass@ \leq 6%
Precision	0.86% Pass@ \leq 10%

Independent test measured absorbance calibration for the unit based on a 3-plate average. Each pipettor is calibrated and reported on. Report included with each unit.

Tools For Life Science Discovery




800.767.0665
Sales@PhenixResearch.com
www.PhenixResearch.com

BenchTop Plate Pipettor

Independent Test Results

Well Volume Goal: 250µl

	1	2	3	4	5	6	7	8	9	10	11	12
A	251.24	251.96	252.06	251.52	251.26	249.54	252.55	251.73	251.37	248.32	248.90	249.40
B	250.26	251.81	251.55	252.40	250.82	252.26	250.99	251.97	251.70	250.50	248.15	248.37
C	248.97	251.34	249.34	248.71	249.81	249.71	250.09	248.76	250.59	249.46	250.75	249.93
D	249.30	249.57	250.99	249.26	249.79	248.62	250.18	246.53	250.21	250.06	248.74	250.75
E	248.05	250.58	250.17	250.41	250.67	249.26	249.92	250.59	249.27	250.31	250.73	246.93
F	250.13	250.66	250.89	250.51	249.19	251.56	249.61	250.08	250.28	249.75	250.72	250.43
G	249.89	250.88	249.57	250.75	248.44	250.83	249.67	248.50	252.07	248.16	250.01	249.47
H	249.10	251.20	250.62	249.61	248.26	249.05	249.79	248.80	251.46	246.55	245.87	249.67



Overall Results: Volume 250µl

Avg. Volume	250.01µl
Standard Deviation	1.31µl
Accuracy	0.00% Pass@≤4%
Precision	0.52% Pass@≤2%

Independent test measured absorbance calibration for the unit based on a 3-plate average. Each pipettor is calibrated and reported on. Report included with each unit.

Tools For Life Science Discovery



800.767.0665
Sales@PhenixResearch.com
www.PhenixResearch.com